

Care
guide

VENACTIF®

Compression Therapy



MADE EXCLUSIVELY IN FRANCE
Feel the difference!



FOLLOW OSSUR ON



WWW.OSSUR.CO.ZA

Ossur South Africa
Units 4 & 5, 3 on London
Brackenridge Business Park
Brackenfell
Cape Town 7560

TEL: 0860 888 123
FAX: 021 981 8925
info@ossur.com




VENACTIF® COMPRESSION HOSIERY FOR A HEALTHIER LIFE

When worn on a daily basis, compression hosiery has proven effective in the treatment of heavy or painful legs. It compresses the veins, facilitating the flow of blood from the feet to the heart and providing instant relief.

For it to be effective, compression hosiery must be worn every day.

WEARING

Putting on your compression hosiery requires a slightly different approach to traditional hosiery.



- A Roll the product down to the heel. Then steadily push your foot in until the heel is firmly in place.
- B Gradually unroll the stocking or sock along the leg by way of light, even movements. For increased comfort, tug the material gently from right to left in order to distribute it evenly and avoid folds.

VENACTIF TIPS

To make fitting easier, we suggest that you thoroughly dry your skin before putting on your compression hosiery.

To prevent tears and ladders (caused by nails, rings, etc.), wear simple plastic gloves, such as household cleaning gloves.

To get rid of folds and creases, try rubbing with gentle, upward movements of your hand.

WASHING

To ensure your VENACTIF® product retains its elasticity and doesn't slip, we recommend washing after each and every use.

Regular washing does not damage the hosiery – it actually prolongs its lifespan.

- Wash on 'delicate wash' cycle at 40° C – wash daily and preferably in a mesh bag
- Use mild detergent and NO fabric softener to protect the fibres and silicone on the anti-slide band
- Dry on a flat surface away from heat sources – no clothes racks, radiators or tumble driers etc
- Do NOT iron or dry clean – both will compromise the product significantly



VENACTIF TIPS

Avoid letting the silicone anti-slide band come into contact with any greasy substance such as moisturising cream or oil.

Wipe the silicone band once a week with cotton wool soaked in alcohol.

